

Antenna Used with QCA6234 in P20830-xx, FujiFilm SonoSite iViz Ultrasound Machine

Antenna: Fractus FR05-S1-NO-1-004

Gain: 1.8 dBi in 2.4-2.5 GHz, 4.9 dBi in 4.9-5.875 GHz

Description: See user manual for device

FCC Part 15.203 compliance: This antenna is unique to this device. It is installed at the manufacturer and is not user serviceable. The design of this is for the antenna to be permanently attached inside of the device and not able to be removed or serviced unless at an authorized service facility.